

# DPM1B

## MEDIDOR DE TRANSDUCTORES NEUMÁTICOS



### Descripción

El DPM1B es un medidor de transductores neumáticos diseñado para medir presiones positivas y negativas de los dispositivos médicos. Tanto en forma de gas como líquida, y para generar presión dentro del rango de  $\pm 300$  mmHg para ayudar en la reparación y en el control de calidad.

### Características principales

- Funcionamiento a pilas
- Genera y mide presiones positivas o negativas
- Funciona con gases o líquidos
- Exactitud dentro del 1%

Tel: 935 862 747  
 Fax: 935 862 749  
 E-mail: [info@stelectromedicina.es](mailto:info@stelectromedicina.es)  
 Web: [www.stelectromedicina.es](http://www.stelectromedicina.es)

**ST - Electromedicina, s.a.**  
 c/ Atenas, nº 23-25, Nave 2  
 Pol. Ind. Cova Solera  
 08191 - RUBI  
 (Barcelona) - SPAIN

## Technical Specifications

### Pressure Measurement

Operating Range:  $\pm 300$  mmHg  
Pressure Generation Range:  $\pm 300$  mmHg  
Accuracy:  $\pm 1\%$  of reading or  $\pm 1$  mmHg  
Units of Measure: mmHg

### Environmental Requirements

Operating Temperature:  $10\text{ }^{\circ}\text{C}$  to  $+40\text{ }^{\circ}\text{C}$   
( $+50\text{ }^{\circ}\text{F}$  to  $+104\text{ }^{\circ}\text{F}$ )

### General Information

Display/Control: 0.5 in LCD with LO BATT indication  
Power: 9 V alkaline battery  
Weight: 10 oz (260 g)

# FLUKE®

## Biomedical

## Ordering Information

### Model

**2249779:** DPM 1B Pneumatic Transducer Tester

### Standard Accessories

**2572314:** User/service manual  
**2242666:** 4-way stop cock adapter  
**614487:** 9 V alkaline battery

### Optional Accessories

**2242666:** 4-way stop cock adapter  
**2212678:** Acrodisc filter

### About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

### Fluke Biomedical Regulatory Commitment

As a medical device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 certified and our products are:

- FDA Compliant
- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required

### Fluke Biomedical.

*Better products. More choices. One company.*

#### Fluke Biomedical

PO Box 9090, Everett, WA 98206-9090 USA

#### Fluke Biomedical Europe AS

Vegamot 8, N-7048 Trondheim, Norway

#### For more information, contact us:

In the U.S.A. (800) 648-7952 or  
Fax (425) 446-5629  
In Europe/M-East/Africa +47 73954700 or  
Fax +47 73954701  
From other countries +1 (425) 347-6100 or  
Fax +1 (425) 446-5629  
Email: [sales@flukebiomedical.com](mailto:sales@flukebiomedical.com)  
Web access: <http://www.flukebiomedical.com>

©2006 Fluke Biomedical. All rights reserved.  
Specifications subject to change without notice.  
Printed in U.S.A.  
LIT047